

Making composite images from layered morphology 2D images

Evelyn Metzger

Vikram Kohli

2024-06-12

Contents

1	Description	1
2	Where to find the script?	2
3	Required libraries	2
4	User inputs (required)	3
5	Additional inputs (changed within the script)	3
6	Output	3
7	Usage	4
8	Example	4

1 Description

Composite images of CosMx™ Spatial Molecular Imager (SMI) fields of view (FOVs) can be useful when using open-sourced software such as [squidpy](#) and [giotto](#). In this post, we describe the `make_composite.py` script, a developmental python script that creates such composite images from layered morphology 2D images that can be exported from the AtoMx™ Spatial Informatics Portal (SIP). Layered images are extracted from the 2D morphology TIF files and written in a file format selected by the user. The extracted images are converted to 8bit, and composite images are written from these 8bit images.



Figure 1: A composite image created using the `make_composite.py` script. This image represents all channels of a single field of view in the [publicly available](#) mouse coronal hemisphere FFPE dataset.

i Note

`make_composite.py` is a development version. Use at your own risk.

2 Where to find the script?

The script and license can be found in the [assets/make-composite](#) folder of the repository.

3 Required libraries

The script requires the following libraries to be installed:

- Pillow (*e.g.*, `pip install pillow`)
- Numpy (*e.g.*, `pip install numpy`)

4 User inputs (required)

- **clipping** (int or float) - Histogram clipping percentage. This value is the percentage of the histogram to clip on the left and right side. The effect changes the contrast of the image. A higher percentage produces more contrast. The user needs to determine the appropriate percentage by testing on a subset of images. The same clipping value is applied to all images. Generally, setting the value between 1 and 3 is a good starting point. Specifying a clipping value of 0 will not alter the histogram. A value is required, there is no default. **clipping** is a required input.
- **user_format** (str) - File format to be written. Options are jpg, png, and tif. All output files will be written in the format specified by the user. **user_format** is a required input.

5 Additional inputs (changed within the script)

- Variable: **colors** = ['cyan', 'red', 'yellow', 'blue', 'magenta']. The variable is the composite color scheme (not a user input; changeable within the script). The colors are listed in order of channel number (channel 0 to channel 4). Example: Channel 0 is colored 'cyan',
- Variable: **compress_value** (set to 3). Lossless file compression value. Higher values produce smaller files at the expense of increased script execution time. The set value is a compromise between file size and execution time.

6 Output

- **raw** - The extracted tif files from the morphology 2D images will be saved in this folder. The file format will follow <fov_num>_ch<#>_raw.<user_format>. <fov_num> is the fov number, <ch<#> is the channel number (from 0 to 4), and <user_format> is the specified file type format (see User inputs). Note: If the specified **user_format** is jpg, the raw files will be 8bit jpg files.

Example

F001_ch0_raw.jpg (for **user_format** = jpg)

- **8bit** - The images in the raw_folder are converted to 8bit and saved in this folder. Note: If the specified **user_format** is jpg, the 8bit files are identical to the raw files. The file format will follow <fov_num>_ch<#>_8bit.<user_format> except when **user_format** = jpg

Example

```
F001_ch0_8bit.tif (for user_format = tif)
F001_ch0_raw.jpg (for user_format = jpg)
```

- **8bit_autocontrast** – Images in the 8bit folder are autocontrasted based on the user supplied clipping value. The file format will follow <fov_num>_ch<#>_8bit_autocontrast.<user_format>

Example

```
F001_ch0_8bit_autocontrast.png (for user_format = png)
```

- **composite** - Composite images created from the images in the 8bit folder. The composite type is a screen composite. The file format will follow <fov_num>_composite.<user_format>

Example

```
F001_composite.jpg (for user_format = jpg)
```

- **composite_autocontrast** - Images in the composite folder are autocontrasted based on the user specified clipping value. The file format will follow <fov_num>_composite_autocontrast.<user_format>

Example

```
F001_composite_autocontrast.png (for user_format = png)
```

7 Usage

```
cd to/your/Morphology2D folder
python /path/to/your/make_composite.py # and follow the on-screen prompts
```

Regex pattern matching on 2D morphology file name format is implemented, however, only NanoString 2D morphology files should be present in the folder containing the make_composite script.

8 Example

The example dataset that we used was the mouse coronal hemisphere FFPE dataset that is available to download from NanoString's website [here](#).

The Morphology2D folder is found within the CellStatsDir folder and has TIF files for each of the 130 FOVs.

```
# In Terminal
cd /path/to/slide/CellStatsDir/Morphology2D
```

```
# In Terminal
```

```
tree -L 1
```

```
20230406_205644_S1_C902_P99_N99_F001.TIF
20230406_205644_S1_C902_P99_N99_F002.TIF
...
20230406_205644_S1_C902_P99_N99_F129.TIF
20230406_205644_S1_C902_P99_N99_F130.TIF
```

Once in the Morphology2D folder, simply run the script and follow the on-screen prompts (Figure 2).

```
# In Terminal
```

```
python /path/to/CosMx-Analysis-Scratch-Space/assets/make-composite/make_composite.py
```

```
*****
Composite script started
*****
Shared Computer: 67%
...
Please specify a clipping percentage as an integer or float: 3
Please specify one of the allowed image format types [jpg, png, tif]: png
Setting the clipping value for autocontrast to 3.0% and the export format to png
Checking for existing folders and creating the necessary folders..
Folders successfully created!
Please wait: Processing the files in /Volumes/Extreme_Pro/data/Helfrain/CellStatsDir/Morphology2D
...Extracting png files (raw) from layered 2D morphology image files
...Converting raw files to 8bit of type png and creating 8bit autocontrast files
/Users/evelymetzger/Documents/Projects/repos/CosMx-Analysis-Scratch-Space/assets/make-composite/make_composite.py:161: RuntimeWarning: invalid value encountered in divide
reduced_bit = (array / array.max())*255
/Users/evelymetzger/Documents/Projects/repos/CosMx-Analysis-Scratch-Space/assets/make-composite/make_composite.py:162: RuntimeWarning: invalid value encountered in cast
img = Image.fromarray(reduced_bit.astype('int8'))
...Writing composite and composite autocontrast images as png
Finished. Total run time = 1558 sec
```

Figure 2: Screenshot of standard output from terminal following script execution. In this example, I set the clipping percentage to 3 and the output to png. On a Macbook Pro M1, this took about 25 minutes to process 130 FOVs.

When complete, the structure of the Morphology2D folder should resemble this:

```
# In Terminal
```

```
tree -L 1
```

```
20230406_205644_S1_C902_P99_N99_F001.TIF
20230406_205644_S1_C902_P99_N99_F002.TIF
...
20230406_205644_S1_C902_P99_N99_F129.TIF
20230406_205644_S1_C902_P99_N99_F130.TIF
8bit
```

```
8bit_autocontrast  
composite  
composite_autocontrast  
raw
```

These composite images can now be imported into open-sourced software or explored further.